SYNMEC

5XFJ-L Laboratory Seed Grading Machine

Separate seeds by seeds width or thickness.









SGS



5XFJ-L Laboratory Seed Grading Machine



Parameter reference:

| Name | Laboratory seed grader |
|--------------------------|------------------------|
| Model | 5FJ-L |
| Capacity | 50 kg/h |
| Motor power | 0.13 kw |
| Rated rotation speed | 1400 r/min |
| Vibration frequency | 0-400 times/min |
| Vibration tank amplitude | 15mm |
| Dimension | 1000×750×1520 mm |

Function:

Laboratory seed grader is used for a small amount of seeds grading.

This quite fits for cereal seeds, grass seeds, flower seeds, vegetable seeds, herbs and seeds and other granular materials separating or dimensions classification.